

Reprint from

RECENT ADVANCES
IN DOPING ANALYSIS
(2)

M. Donike
H. Geyer
A. Gotzmann
U. Mareck-Engelke
(Editors)

Sport und Buch Strauß, Köln, 1995

C. Ayotte, A. Charlebois, D. Goudreault, C. Lévesque:
Improvements of Drug Detection and Analysis by Appropriate Preparation of the Injection
Port of the Gas Chromatograph
In: M. Donike, H. Geyer, A. Gotzmann, U. Mareck-Engelke (eds.) Recent advances in doping
analysis (2). Sport und Buch Strauß, Köln, (1995) 375-392

C. Ayotte, A. Charlebois, D. Goudreault, C. Lévesque.

Improvements of drug detection and analysis by appropriate preparation of the injection port of the gas chromatograph.

Institut national de la recherche scientifique (INRS-Santé), Montréal, Québec.

In this paper, we will present several examples of drastic improvements of the quality of results obtained in GC/MS by reducing the "reactivity" of the injection chamber. The carrier gas (zero grade) is purified (O_2, H_2O) with a OMI-1 Indicating Purifier (Supelco, Canada) containing a resin of solid lithium compound. The glass liner and the glass wool are silanized overnight in a 5% DMDCS/toluene solution and rinsed three times with methanol. One of the most striking result was the obtention of the TMS-derivative of pemoline¹ described as impossible to detect through chromatography.

In the following pages, we are presenting selected examples taken from the GC/MS analysis of some TMS-ether, TMS-enol derivatives of anabolic androgenic steroids or their metabolites : oxandrolone, 6 β -hydroxymethandienone, 6 β -hydroxydehydrochlormethyltestosterone, stanozolol, T/E values and epi trenbolone. The improvements are observed with the detection levels, long term stability and high reproductibility.

The examples of propoxyphene, strychnine and codeine are also given.

¹ C. Ayotte: Old dope agent, old method, 11th Köln Workshop in dope analysis, Köln (1993).

Procedure IV

GC/MS analyses in the SIM (selected ion monitoring) mode are carried out with Hewlett Packard HP 5890 gas chromatographs with direct coupling to HP-MSD 5970 quadrupolar filters. Systems are equipped with automatic samplers 7673 A and controlled by HP 59970 and Unix MS Chemstations. The injections are carried out in the splitless mode (1 μ L) and the separation is achieved on HP-5 capillary columns.

Chromatographic parameters

carrier gas : He

injection port : 270°C

transfer line : 310°C

injection mode : splitless 30 sec.

initial temperature of the oven : 100°C (1 min)

initial rate : 16°C

first temperature : 220°C

final rate : 3.8°C

final temperature : 300°C (10 min)

PROCEDURE IV

SCREENING/CONFIRMATION

(Free and glucuronides fraction)

2 to 10 mL urine

↓
add

5 µg epietiocholanolone
500 ng 17α-methyl-5α-androstan-3β,17β-diol

↓

apply on Sep-Pak C18 cartridge
wash with 5 mL of water
wash with 3 mL of n-hexane
elute with 5 mL of methanol
evaporate methanol to dryness

↓

add

1 mL phosphate buffer 0.2M pH: 6.9
β-glucuronidase from E. Coli (type IX-A, 2000 units)

↓

37°C for 16 hours or 50°C for 1 hour

↓

add

100 mg solid buffer pH 9
5 mL diethyl ether

↓

shake mechanically 10 min. and centrifuge 5 min.
evaporate organic phase to dryness

↓

dissolve residue in 150 µL of methanol
transfer to vial containing 500 ng of 5α-androstan-17-one

↓

evaporate to dryness

↓

add


2µL ethanethiol
100 µL MSTFA
2 µL TMSI

↓


70°C for 30 min.

↓


Inject into GC/MS and GC/FID




GC/MS (SIM mode) analysis of the TMS-ether, TMS-enol derivatives of a mixture of standards. Reactives used for the chemical derivatization; MSTFA:TMSI: ethanthiol (upper panel), MSTFA:TMSI (lower)




GC/MS (SIM mode) analysis of a mixture of standards injected over a 18 hours period on a silanized glass liner.




GC/MS (SIM mode) analysis of a mixture of standards injected over time
(enlarged)



GC/MS (SIM mode) analysis of a mixture of standards (first, third and tenth injections)



GC/MS (SIM mode) analysis of a mixture of standards (first, third and tenth injections) (enlarged)




Calibration Table


[T]/[Et] Standards spikes/urine (Solutions #080394)
 Calibration Table Last Updated: Thu Apr 07 14:07:04 1994
 Reference Peak Window: 5.000 % of Retention Time
 Non-reference Peak Window: 5.000 % of Retention Time
 Default Sample Amount: 0
 Uncalib. Peak Response Factor: 0
 Default Multiplier: 1

Ret Time	PK#	Signal Descr	Amt ng/ul	Lvl	Area	Pk-Type	Partial Name				
23.465	1	432.40 amu	0.05007000	3	48448		Epitestostero				
			0.1001000	1	92560						
			0.2003000	4	199001						
			0.3004000	2	251919						
			1.001000	7	867146						
			3.004000	6	2550640						
			5.007000	8	4514220						
			10.01000	5	5568990						
			24.425	2	432.40 amu	0.05043000		1	66190		Testosterone
						0.1009000		2	106023		
0.5043000	3	403990									
1.009000	5	625639									
0.8069000	4	736529									
3.026000	6	2551990									
6.052000	7	4944800									
10.09000	8	8955800									
25.173	3	435.40 amu				5.000000	5	903618	+I	Istd 2	
						5.000000	3	1068720			
			5.000000	1	1150270						
			5.000000	2	1274890						
			5.000000	4	1278840						
			5.000000	6	1315420						
			5.000000	7	1315760						
			5.000000	8	1425250						


Calibration curves used for the determination of both testosterone and epitestosterone concentrations. 17 α -methylandrostan-5 α -3 α ,17-diol is used as internal standard. Analytical standards are spiked in urine samples and



Ion 432.4 and 435.4 amu from the GC/MS analysis of the control urine sample injected over a 18 hours period each 4 samples.





Superimposition of three ion chromatograms (432.4, 435.4 amu) obtained from the analysis of the control urine sample injected over a 18 hours period.





Superimposed ion chromatograms (ions 412.4 and 307.3) obtained from the analysis of the control urine sample over a 18 hours time period.

File : C:\HPCHEM\1\DATA\CAP92\ADT30BL.D
 Operator : d.g.
 Acquired : 15 Jan 92 3:20 pm using AcqMethod BLEUE.M
 Instrument : MS_5971
 Sample Name: adt-30 proc 1derive 20ul
 Misc Info :
 Vial Number: 1









GC/MS (full scan mode) analysis of a propoxyphene urine sample
 injected on an untreated glass liner

File : C:\HPCHEM\1\DATA\SMCC\ADT-30.D
Operator : C.L.
Acquired : 6 Mar 94 8:46 pm using AcqMethod STIM2
Instrument : 5971 - In
Sample Name: ADT-30 CAP 92 PR.I
Misc Info : PROPOXYPHENE
Vial Number: 8



GC/MS (full scan mode) analysis of a propoxyphene urine sample injected on a silanized glass liner.

Fichier: C:\HPCHEM\1\DATA\SMCC\STD1.D
 Opérateur: MANON
 Date analysé: 6 Mar 94 12:34 pm
 Instrument: 5971 - In
 Méthode utilisée: SIM.M
 Échantillon: STD DE PERFORMANCE
 Information suppl: WITHOUT SILANIZED GLASS WOOL
 Vial no: 2









RÉSULTATS

Produit	RT	Surface	Ratio
MA	8.99	57183377	3.07
PDM	11.95	43919861	2.36
MPH	14.41	38487621	2.07
Peth	14.57	18625288	
Codéine	22.38	3452812	185.38

GC/MS (SIM mode) analysis of a mixture of standards injected on silanized glass liner.

Fichier: C:\HPCHEM\1\DATA\SMCC\STD.D
 Opérateur: MANON
 Date analysé: 6 Mar 94 10:48 am
 Instrument: 5971 - In
 Méthode utilisée: SIM.M
 Échantillon: STD DE PERFORMANCE
 Information suppl: SILANIZED GLASS WOOL
 Vial no: 2





RÉSULTATS

Produit	RT	Surface	Ratio
MA	9.00	51597207	3.06
PDM	11.95	36445071	2.16
MPH	14.41	40323176	2.39
Peth	14.56	16871694	
Codéine	22.36	7242364	429.26



GC/MS (SIM mode) analysis of a mixture of standards injected on silanized glass wool and liner.

File : C:\HPCHEM\1\DATA\SMCC\STRYC1.D
 Operator : MANON
 Acquired : 6 Mar 94 12:04 pm using AcqMethod STRYC
 Instrument : 5971 - In
 Sample Name: STRYCHNINE 0.5UG/ML
 Misc Info : WITHOUT SILANIZED GLASS WOOL
 Vial Number: 1



GC/MS (full scan mode) analysis of strychnine authentic standard injected on silanized glass liner.

File : C:\HPCHEM\1\DATA\SMCC\STRYC.D
 Operator : MANON
 Acquired : 6 Mar 94 10:17 am using AcqMethod STRYC
 Instrument : 5971 - In
 Sample Name: STRYCHNINE 0.5UG/ML
 Misc Info : LINER WITH SILANIZED GLASS WOOL
 Vial Number: 1



GC/MS (full scan mode) analysis of a strychnine authentic standard injected on silanized glass wool and liner.